



BEDS

Accessories

clinic



clinic design

www.bettenmalsch.de




Standard handle




Handle with automatic winder



Intravenous drip clasp for trapeze bars




Model	Bottle holder
straight	plastic/stainless steel
curved	plastic/stainless steel



Mattress wedge for integrated bed extension



5th castor




Additional medical rail




Urine bottle and basket



Rescue blanket



Mediline 5 zone mattress with edge reinforcement



Viscoline 2 layer mattress with edge reinforcement



Oxygen bottle holder



Flexible hand controller holder



Second patient hand control (cross) with clip



Bed transport trolley



Bedding storage




Under-bed lighting



Name sign holder



Restraint holder



Hyper Foam anti-decubitus mattress

"WXT" SERVER (including static tray)

Static tray 18 cm x 40 cm



Magazine rack optional available

"WX" SERVER

Tray width: 60 cm

Tray depth: 40 cm

Sub frame width: 72 cm

Sub frame depth: 42 cm

Height adjustment 72 cm – 111 cm

Special features

- ▶ Moving frame made from precision tubing, powder coated
- ▶ 4 furniture castors Ø 5 cm, of which 2 are lockable
- ▶ Aluminium column with cushioned raising / lowering mechanism
- ▶ Laminate tray with lip on both sides, variable, both sides can be tilted



"WK I"

BEDSIDE CABINET

Width: 45 cm

Depth: 49 cm

Height: 55 cm

"WK I-N"

BEDSIDE CABINET

Width: 45 cm

Depth: 49 cm

Height: 77 cm

Arrangement:

- 1 drawer – lockable as an option
- 1 deep drawer with bolder holder or revolving door, with shelf
- 1 open storage compartment (WK I-N)
- 4 furniture castors Ø 7.5 cm, of which 2 are lockable



"WK III" BEDSIDE CABINET

Width: 45 cm

Depth: 49 cm

Height: 77 cm

Tray height adjustment: 77 cm - 115 cm

Arrangement can be used on both sides:

- 1 drawer with push-through limiter
- 1 side folding bed tray 60 cm x 40 cm
- 1 open storage tray
- 1 deep drawer with bottle holder or revolving door, with shelf

Special features

- ▶ 5 furniture castors Ø 7.5 cm, of which 2 are lockable
- ▶ Aluminium column with cushioned raising / lowering mechanism
- ▶ Laminate tray with lip on both sides, variable, both sides can be tilted

